# **F-SAFE Sicherheitsschränke / F-SAFE Safety Cabinets / F-SAFE Armoires de sécurité**

## **Montage- und Bedienungsanleitung / Installation and Operating Manual / Instructions de montage et d'utilisation**

SST-P 12/20 F30  
SST-P 12/20 F90

SST-P 6/20 F30  
SST-P 6/20 F90

SST-P 11/6 F30  
SST-P 11/6 F90

SST-LI 12/20 F90

SST-LI 6/20 F90

SST-LI 11/6 F90

Version: 02/19  
EP.V.9021

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Dear customer,

Thank you for choosing a premium PROTECTO product!

Please read this operating manual carefully to ensure that you always get the best out of the system as a whole. Failure to observe the information provided here will render the warranty null and void.

We expressly reserve the right to make technical changes aimed at improving the housings, or which are due to legal changes, without prior notice.

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Reprinting or reproduction of this Installation and Operating Manual, or of parts thereof, is subject to the prior written permission of PROTECTO.

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**DEUTSCH**

## Contents

<b>1. General safety instructions</b>	38
<b>2. Technical data</b>	
2.1. Sales and customer service	39
2.2. Description / use	39
2.3. Type of construction	39
2.4. Technical data	40
<b>3. Transport</b>	
3.1. Safety precautions for transport	41
<b>4. Installation</b>	42
<b>5. Commissioning</b>	
5.1. Wing doors	43
5.2. Installation of the shelves	43
5.3. Sump / perforated plate insert	44
5.4. Replacing / checking the sump	44
5.5. Storage volumes	44
<b>6. Ventilation</b>	
6.1. General	44
6.2. Connection to an exhaust system	45
6.3. Selection of the technical ventilation	45
6.4. Operation without technical ventilation	46
<b>7. Operation</b>	
7.1. General principles	46
7.2. Instructions on storing lithium ion batteries (SST-LI)	46
<b>8. Decommissioning</b>	47
8.1. After a fire	47
8.2. Disposal	47
<b>9. Maintenance</b>	48
9.1. Daily inspection	48
9.2. Monthly inspection	48
9.3. Inspection by our maintenance service	48
9.4. Troubleshooting	49
<b>10. Operating information</b>	49
<b>11. Connection diagrams</b>	50
<b>12. Optional accessories</b>	
12.1 General safety information	52
12.2 Intended use	53
12.3 Warranty, surrender of goods	53
12.4 Single-cable / single-pipe penetration	54
12.5 Cable / pipe bulkhead	56
12.6 Recirculation fan, type SST-UL V4	59
12.7 Exhaust fan, type SST AL	65
<b>13. General construction approval / GS mark certificate (German)</b>	105

## 1. General safety instructions

- Observe the laws and regulations applicable to the handling of hazardous substances and the instructions in this operating manual.
- The operating/ assembly instructions must always be kept at the place of operation
- Observe all accident prevention regulations and the Workplace Regulation.
- Legally required safety instructions are to be posted in a visible place for users.
- All safety and hazard warnings and the rating plate must be kept in a readable condition.
- Prior to commissioning, a functional test must be carried out to detect possible damage due to transport, shipment or improper installation.
- Do not use the product if it is not in proper condition.
- Observe the technical data (storage quantities, loads, etc.) for the respective product.
- Improper use may impair the protective function.
- Operating personnel must be instructed in advance to prohibit access by unauthorised persons.
- The required safety checks may only be carried out by authorised specialist personnel.
- Always keep the opening areas of the doors free.
- Doors / drawers should only be kept open for loading and unloading.
- Ensure sufficient technical ventilation (see Section 6).
- Check the resistance of cabinet surfaces before storing chemicals.
- Absorb and dispose of escaping hazardous substances immediately.
- For the storage of aggressive substances (acids, alkalis, corrosive gases), special acid and alkali cabinets should be used – risk of corrosion of interior fittings and ventilation components.
- It is beyond the scope of this document to list all applicable regulations and guidelines. Always verify that your system complies with applicable regulations.

## 2. Technical data

### 2.1. Sales and customer service

Protectoplus Lager- und Umwelttechnik GmbH  
Grüner Kamp 19-21  
24768 Rendsburg  
Germany  
Internet: [www.protecto.de](http://www.protecto.de)  
Email: [info@protecto.de](mailto:info@protecto.de)

### 2.2. Description / use

Safety cabinet according to DIN EN 14470, Part 1 (date: 2004-04)  
PROTECTO safety cabinets are used to store flammable liquids in workrooms in conformity with applicable national regulations.

### 2.3. Type of construction

#### 2.3.1. Body:

- 2.3.1.1. Type F30 fire protection panels
- 2.3.1.2. Type F90 fire protection panels with surface coating, material classification A2, non-combustible

#### 2.3.2. Door closer:

##### 2.3.2.1. Door closer 1:

- Automatic door closing by means of an oil-damped gear drive closer
- Temperature-dependent central locking device as standard
- Included with:  
SST-P 11/6 F30, SST-P 6/20 F30, SST-P 12/20 F30,  
SST-P 6/20 F90 (until 2010), SST-P 12/20 F90 (until 2010)

##### 2.3.2.2. Door closer 2:

- Automatic door closing in case of fire
- Included with:  
SST-P 11/6 F90, SST-P 6/20 F90, SST-P 12/20 F90,  
SST-LI 6/20 F90, SST-LI 12/20 F90, SST-LI 11/6 F90

#### 2.3.3. Features:

Height-adjustable shelves, one (1) removable perforated plate insert as the first storage level, one sump

## 2.4. Technical data

<b>Type F30 Fire resistance 30 minutes</b>				
Model		EN32.196.120	EN31.196.060	EN32U.063.110
Type		<b>SST-P 12/20 F30</b>	<b>SST-P 6/20 F30</b>	<b>SST-P 11/6 F30</b>
External dimensions				
Height	approx. mm	1960	1960	635
Width	approx. mm	1195	595	1100
Depth	approx. mm	595	595	550
Empty weight with shelves and sump	approx. kg	226	147	97
Max. total weight when fully loaded	approx. kg	526	297	172
Usable volume of safety cabinet	approx. m <sup>3</sup>	0.88	0.39	0.17
Capacity of sump	approx. l	44	19	24
Capacity of largest individual container	approx. l	35	15	18
Shelf loading capacity for evenly distributed loads	approx. kg	75	75	75

<b>Type F90 Fire resistance 90 minutes</b>				
Model		EN92.196.120	EN91.196.060	EN92U.063.110
Type		<b>SST-P 12/20 F90</b> <b>SST-LI 12/20 F90</b>	<b>SST-P 6/20 F90</b> <b>SST-LI 6/20 F90</b>	<b>SST-P 11/6 F90</b> <b>SST-LI 11/6 F90</b>
External dimensions				
Height	approx. mm	1960	1960	670
Width	approx. mm	1200	600	1100
Depth	approx. mm	595	595	500
Empty weight with shelves and sump	approx. kg	600	340	210
Max. total weight when fully loaded	approx. kg	900	640	290
Usable volume of safety cabinet	approx. m <sup>3</sup>	0.83	0.38	0.18
Capacity of sump	approx. l	44	20	24
Capacity of largest individual container	approx. l	40	18	21
Shelf loading capacity for evenly distributed loads	approx. kg	75	75	75

## 3. Transport

### 3.1. Safety precautions for transport

Bear in mind that the cabinets are heavy; a minimum of 2 persons are required to transport them.

**Important:**

- **The transport locks in the door joints must not be removed until the cabinets reach the place of use!**
- **Type SST 12/20 F90 and SST 11/6 F90 safety cabinets as well as type SST-LI 12/20 F90 and SST-LI 11/6 F90:  
The transport lock in the rear wall must not be removed until the cabinets reach the place of use!  
You must remove this lock to ensure that the doors are locked in case of fire.**

The PROTECTO safety cabinet may only be lifted for transport using suitable means, which should be attached to the underside; it must always be transported upright. Be careful not to tilt it during lifting or transport. Tilting could damage the cabinet and impair or destroy the functionality of the built-in safety features. The cabinet should be set down again gently at the place of use. All floor standing cabinets are shipped with a base and can be lifted by the door after removing the base panels.


The 600 mm wide floor standing cabinet can only be lifted using an industrial truck with a fork. There is a risk of toppling over if the cabinet is handled incorrectly. The base panels should be secured to the cabinet again after reaching the place of use because they are essential for fire protection. Without them, the cabinet's fire integrity will be compromised and the approval will cease to be valid.

## 4. Installation

The PROTECTO safety cabinet must be installed on a level surface which is capable of supporting it when fully loaded. Shims should be placed underneath the cabinet if the floor is uneven. The cabinet should be securely screwed to the building or structure for safety reasons to prevent it from toppling over.

Remove all transport protections. Only then are the products fully ready to operate. Remove the key and the bow handle(s) from the top of the cabinet and fasten them to the door(s) using the threaded screws which are enclosed. The PROTECTO safety cabinet must be absolutely level. You achieve this by adjusting the feet in the base of the floor standing cabinet (they can be adjusted by 25 mm). Insert a hex key (5 mm) through the holes in the bottom of the cabinet from the inside to adjust the feet to the required height. The self-closing function will be impaired if the cabinet is not level.

After aligning, mount the 42mm (EN9x.xxx) or 38mm (EN30.xxx) plinth on the front of the floor standing cabinet using the enclosed 70mm long screws with covers.

After aligning the floor standing cabinet, secure the base panel using the screws provided.

The base panels must be fitted again because they are essential for fire protection. Without them, the cabinet's fire integrity will be compromised and the approval will cease to be valid.

The actual place of use / operation must be such that the safety cabinet

- Cannot be damaged by vehicles,
- Is not exposed to direct or indirect heat,
- Is protected against moisture,
- Is not exposed to operating temperatures less than 5°C or greater than +40°C.

### Instructions on storing lithium ion batteries (SST-LI):

The combustion of Li-ion batteries can lead to heavy smoke. The cabinets primarily serve to contain the fire, but are not smoke-tight. PROTECTO recommends installing a smoke detector close to the cabinet to control on-site ventilation to extract flue gases in case of fire.

## 5. Commissioning

### 5.1. Wing doors

The doors can be locked by means of the integral locking cylinder and opened independently of one another.

The door closer is set at the factory to operate correctly at room temperature. If extreme temperature changes occur in your application, it could take up to 12 hours to adapt to room temperature.

**Important: The closing force and / or speed of the door closer may only be altered by authorised persons.**

The central locking device is released at a temperature of 50°C (+0/-10°C) base and locks the door automatically.

### 5.2. Installation of the shelves

If the cabinet includes tray shelves as special shelves (for discharging leaking liquids into the sump):

- The shelves which are shipped with the cabinet and the shelf support pins (4 or 5 per shelf) and fixing screws (4 per shelf) are installed in the hazardous material cabinet at the factory.
- To position the shelves, attach the shelf support pins to the side walls in the required position.
- Place the shelf on the shelf support and screw the fastening screws (4 per shelf) through the tray shelf into the side wall of the cabinet (tilt protection and safety device).

Important for equipotential bonding: when fitting the screws at the rear right, the locking washer must be placed between the exterior of the special shelf and the silver foil.

If the cabinet includes a base pan as a special shelf (for collecting leaking liquids):

- Screw the base pan to the side walls using four 4.5 x 35 mm screws to prevent it from tilting or being pulled out inadvertently.

The top shelf must not be more than 1.75 m above the floor.

### 5.3. Sump / perforated plate insert

The sump in the floor of the PROTECTO safety cabinet serves to collect leaking liquids inside the cabinet; it should not be used as additional storage space for containers.

If the sump is fitted with a removable perforated plate insert, on the other hand, the plate may be used for this purpose. In this case, the capacity of the sump is reduced by about 2%.

If the cabinet includes a base pan as a special shelf, the retention capacity specified under "Technical data -> Capacity of the sump" applies.

### 5.4. Replacing / checking the sump

Any leaked liquids which have collected in the sump should be removed and properly disposed of prior to removing the pan. Remember to connect the sump to the potential connection again after removing it.

### 5.5. Storage volumes

The sump must have a minimum retention capacity equivalent to 10% of the (total) storage volume or at least 110% of the capacity of the largest individual container. Please refer to the marking on the outside of the cabinet. This depends on the circumstances of your application: you may need to install a bigger sump if the capacity of the largest individual container is greater than that of the sump for your particular cabinet type.

## 6. Ventilation

### 6.1. General

All models can be connected to an technical exhaust system via the exhaust air outlet. The required air change intervals, maximum pressure drop and exhaust air ducts depend on the applicable national regulations.

Example: Germany

According to TRbF 20, Annex L (Technical Regulations for Flammable Liquids), PROTECTO safety cabinets should, for safety reasons, be connected to a technical exhaust system which leads into the open at a place where no danger exists. At least 10 air changes per hour must take place inside the ventilated safety cabinet when the doors are closed, and the pressure drop must not exceed 1.5 mbar (150 Pa). The air changes should be checked and documented, e.g. using an air flow test tube!

If the PROTECTO safety cabinet is connected to an existing technical ventilation system, the ventilation must operate for 24 hours a day.

Model	Type	Flow rate 10 air changes	Pressure drop in cabinet
EN31.196.060	SST 6/20 F30	approx. 5.0 m <sup>3</sup> /h	approx. 1 Pa
EN32.196.120	SST 12/20 F30	approx. 10.0 m <sup>3</sup> /h	approx. 1 Pa
EN32U.063.110	SST 11/6 F30	approx. 5.0 m <sup>3</sup> /h	approx. 1 Pa
EN91.196.060	SST 6/20 F90	approx. 5.0 m <sup>3</sup> /h	approx. 1 Pa
EN92.196.120	SST 12/20 F90	approx. 10.0 m <sup>3</sup> /h	approx. 1 Pa
EN92U.067.110	SST 11/6 F90	approx. 5.0 m <sup>3</sup> /h	approx. 1 Pa

## 6.2. Connection to an exhaust system

Supply and exhaust air vents are provided on the top and rear of all PROTECTO safety cabinets. The connection fitting should be attached to the exhaust air outlet (see diagram) and the exhaust air duct (metal or flame-retardant plastic pipe) connected to it. Room air can be used as the supply air; under normal conditions of operation no additional, external air supply is required.

The installation of a technical ventilation system, or the connection to an existing exhaust system, does not form part of the PROTECTO portfolio. Please ask a ventilation specialist near you for assistance.

## 6.3. Selection of the technical ventilation (see also section 12.6)

The design and operation of technical ventilation systems which are used in explosive atmosphere are currently subject to European Directive 2014/34/EU (ATEX). Only technical ventilation systems which conform to the technical specifications set out in this directive may be used. The equipment used by us (SST-UL and AL) is designed, manufactured, marked, tested and certified in accordance with the above standard.

No specific hazards for people or the environment arise provided this equipment is used as intended and only with original accessories, and provided this Installation and Operating Manual as well as the information on safe commissioning are observed.

All forms of use which deviate from or exceed the limits of use described above are considered to be contrary to the intended purpose. The manufacturer / supplier is not liable for any damages resulting from such use; the sole responsibility lies with the user.

Before installing the technical ventilation system, check whether your local in-house grid is suitable for operation. The technical ventilation system must be connected and measured by a qualified electrician. The results of the measurements should be documented.

## 6.4. Operation without technical ventilation

The PROTECTO hazardous material cabinet can also be installed and operated in accordance with TRGS 510, Annex 3, Section 2.2 without being connected to a technical ventilation system.

In this case, please note:

- To avoid electrostatic charging, the cabinet must be connected to earth. To do this, connect the standard earth connection on the back of the cabinet to a suitable earthing cable
- The inside of the cabinet is classed as a Zone 1 hazardous area
- The area within a radius of at least 2.5 m around the cabinet up to a height of at least 0.5 m above the floor is a Zone 2 hazardous area
- If the place of use (workroom) has technical ventilation with a minimum of 5 air changes, the Zone 2 hazardous area is reduced to 1 m around the cabinet and 0.3 m above the floor
- The hazardous area should be marked. In addition, all other measures related to an Ex zone must be observed (among other things, no ignition sources, use of Ex-protected devices).

## 7. Operation

### 7.1. General principles

#### Markings

All models are marked according to EN 14470-1 requirements.

#### Storage volumes / joint storage

Observe the applicable national regulations depending on the type classification for your cabinet.

### 7.2. Instructions on storing lithium ion batteries (SST-LI)

When storing lithium batteries, please refer to GDV publication "VdS guideline 3103 (2016-05) - leaflet on damage prevention."

#### Door closure

When storing lithium-ion batteries, the cabinet door must always be kept closed and may only be kept open during loading or unloading.

## 8. Decommissioning

The following time limits should be observed owing to the risk of fire and explosion:

### 8.1. After a fire

Caution: After a fire, the cabinets must not be opened until a period of at least 6 times the duration of the fire has elapsed. Depending on the duration of the fire, an explosive steam-air mixture may have developed. Before opening the cabinet, remove all possible sources of ignition within a radius of at least 10 m. Open the doors extremely carefully. Use only non-sparking tools for this purpose.

Opening after a fire should only be done by the fire brigade or an authorised person.

Provide suitable extinguishing agents for the opening process and wear personal protective equipment.

One indication of whether the waiting time before opening the cabinet must be extended is the temperature of the cabinet surface. If it is still warm, do not open it yet.

After a fire, the PROTECTO safety cabinet may no longer be used. In the event of a fire, it can never be ruled out that the effects of heat may have damaged the fire protection features and that the cabinet has lost its tested fire resistance time of 90 minutes.

### 8.2. Disposal

The cabinet can be sorted for recycling and need not be disposed of as hazardous waste due to its use of harmless materials.

## 9. Maintenance

The PROTECTO safety cabinet is regarded as a safety system. According to Section 4 (3) of the German Workplace Regulations (ArbStättV), Section 10 of the Industrial Safety Regulations and BGR 234 (Rules of the Employers' Liability Insurance Association), a safety audit must be carried out at least once a year by an authorised expert.

The PROTECTO safety cabinet should be inspected for visible damage or defects after installation and prior to commissioning as well as subsequently at the intervals specified below and following any modifications or repairs.

### 9.1. Daily inspection

- Any leaked liquids should be immediately collected and properly disposed of.
- Check your PROTECTO safety cabinet regularly for mechanical damage, fouling and corrosion.

### 9.2. Monthly inspection

The monthly inspection comprises the following:

**9.2.1.** Door closer: Open the door and check that it closes correctly.

**9.2.2.** Ventilation: Check that the ventilation system is working correctly using a wool thread, measuring instrument, etc.

**9.2.3.** Sealing strips: Make sure the sealing strips are correctly seated in the door rebate and front face(s). Replace the sealing strips immediately in case of visible damage.

### 9.3. Inspection by our maintenance service

The inspection of the safety cabinets is laid down in DIN EN 14470, Part 1. Yearly maintenance inspections are the norm.

We offer both single inspections and maintenance agreements.

For more information, please get in touch with us directly on:

Phone: +49 4331/4516-27

Email: [wartung@protecto.de](mailto:wartung@protecto.de)

## 9.4. Troubleshooting

Fault	Possible cause	Possible solution
Door does not close	Safety cabinet out of plumb	Position the safety cabinet so that it is level (by adjusting the feet)
No extraction	Fire protection valve closed because fuse defective	Have the fuse replaced by a customer service technician
Doors stiff	Air exchange rate in cabinet too high, moving parts corroded (hinges, etc.)	Reduce the air change rate, remove rust, oil parts, don't store aggressive substances in the hazardous material cabinet, notify customer service

Please indicate the cabinet type, serial number and model code when reporting any malfunctions to our customer service technicians.

A detailed description of the malfunction will help us identify what caused it more easily.

## 10. Operating information


The owner must compile readily understandable operating information in the language spoken by the employees, taking account of these operating instructions and the specific conditions on site.

The operating information must be posted by the owner at a suitable location at the place of work.


Annual audits of the cabinet, and in particular the heat-triggered release of the central locking device as well as the sealing strips, should be specified in the operating information as a minimum requirement.

## 11. Connection diagrams

### Earthing


### Ventilation



## Ventilation


Type EN31 top view


Type EN91 top view


Type EN32 top view

Type EN92 top view

Type EN32U rear view

Type EN92U rear view

Technical changes are reserved!

## 12. Optional accessories

### 12.1. General safety information

- Keep this Installation and Operating Manual adjacent to the cable / pipe bulkhead for each product.
- Affix any legally prescribed safety information at a location where it is conspicuous to users.
- Strictly observe:
  - All relevant VDE regulations
  - Fire authority regulations
  - General planning and building laws and regulations
  - Accident prevention and workplace regulations
  - Any instructions given by the technical supervisory service
- Pay attention to DIN 4102 regarding the assessment of fire protection and structural aspects.
- Work on electrical systems must be carried out with the power disconnected and performed by electricians.
- Never use a product that is not in perfect condition.
- Ensure that all required safety audits are carried out by an authorised expert.
- Damage or malfunctions due to incorrect transport or installation can only be reliably avoided if these activities are carried out by persons trained and authorised by us.
- Heed the technical information given in our catalogue and data sheet regarding the product concerned.
- Note that “hot spots” may form depending on the position of the fitted components. These components should therefore be distributed as evenly as possible. Check this if necessary by taking measurements inside the cabinet.
- Unauthorised persons must be denied access to the cable / pipe bulkhead.
- Incorrect installation could impair the protective function.
- Observe all information contained in this Installation and Operating Manual.
- Use only suitable, approved screws / dowels to fasten the cabinet to the building or structure.
- Make sure all screws / dowels are suitable for your particular building / structure.
- The doubler plates of the cable / pipe bulkhead must not be glued, painted or otherwise coated.
- There is not enough space in this manual to list every relevant regulation or directive. Please therefore make sure your installation is in conformity with all applicable regulations.

## 12.2. Intended use

- The product described in this Installation and Operating Manual was manufactured in accordance with the latest state of the art as well as generally recognised safety standards
- It may only be used:
  - For its intended purpose, and
  - If in perfect order and condition
- The doubler plates are comprised of coated, non-metallic, plate-shaped components, which are painted light grey (similar to RAL 7035) on both the inside and outside. The edges are in a different colour
- The fire insulation properties are achieved by processing fire protection panels to DIN 4102
- Provided they are installed in the proper manner, the cable / pipe bulkhead has a fire resistance rating of more than 30 or 90 minutes and is smoke-tight in the presence of (smouldering) cable fires
- The information contained in this manual should be strictly observed to ensure that the product and the cable entry are handled correctly
- Max. cable occupancy: 60%
- All other forms of use are considered to be contrary to the intended purpose. PROTECTO is not liable for any direct or consequential damages resulting from such use

## 12.3. Warranty, surrender of goods

All installation instructions must be observed in order to guarantee optimal functioning of our products.

PROTECTOPLUS cannot be held responsible for damages to other products which result from the installation or use of this product.

The warranty which is given by PROTECTOPLUS is limited to products shipped by us.

Design modifications carried out without prior consultation with PROTECTOPLUS result in voiding of the approval / warranty.

Installation services must be guaranteed by the entity responsible for the installation.

Warranty and liability claims for personal injury or property damage are excluded if they are the result of one or more of the following:

- Improper installation, commissioning, operation or maintenance,
- Failure to observe the information on transport, storage, operation and installation,
- Inexpert repairs, or
- Emergencies due to third party actions or force majeure.

### **Industrial property rights**



Utility model documents have been filed with the German patent office to protect the innovation and design.

## 12.4. Installation and operating instructions for single-cable / single-pipe penetration

### 12.4.1. Scope of supply

### 12.4.2. Positioning (if no cable / pipe bulkhead installed)


### 12.4.3. Installation procedure

Data subject to change!

## 12.5. Installation and operating instructions for fire-resistant cable / pipe bulkhead

### 12.5.1. Scope of supply

- 3x doubler plate 200 x 200 mm, with central through-hole Ø135 mm and countersunk fixing holes in the four corners
- 12x fixing / countersunk head screw 4.5 x 70 mm
- 2x firestop plug Ø 134 mm
- 1x firestop silicone cartridge
- 1x Installation and Operating Manual

### 12.5.2. Installation procedure

- The position of the cable / pipe bulkhead (on our products) should be determined on site, taking the dimensions of the doubler plate into account
- Drill a  $\varnothing$  135 mm hole with a core drill (see Figure 1)
- The doubler plates can be used as a marking template

Fig. 1

- Smooth the drilled hole with sandpaper (CAUTION: The edges are sharp)
- Mark and drill  $\varnothing$  3.5 x 35 mm pilot holes on site for fixing the doubler plates (see Figure 2)

Fig. 2

- Position the doubler plates exactly over the  $\varnothing$  135 mm hole, then fasten them using the enclosed countersunk head screws after drilling the pilot holes in the wall
- The fixing holes in the doubler plate are offset, so that one plate can be screwed to another if the plates are turned  $180^\circ$

- After fastening the doubler plates, insert the firestop plugs (one each on the inside and outside (see Figure 3)).

Fig. 3

- Feed the cables through, then seal the openings with the silicone which is supplied with the product.

- F30 cable / pipe bulkheads: Two doubler plates plus the wall plate on which the bulkhead is installed -> wall thickness 120 mm
- F90 cable / pipe bulkheads: Three doubler plates plus the wall plate on which the bulkhead is installed -> wall thickness 150 mm
- Note that the required wall thickness also includes a layer of internal insulation (approx. 40 mm)
- 2 doubler plates: Max. cable diameter which can be fed through = 22 mm
- 3 doubler plates: No max. cable diameter

## 12.6. Recirculation fan, type SST-UL V4

### 12.6.1. Instructions for safe use

To ensure the safe operation of the system, the following basic requirements must be met:

- Temperature at place of installation: between +10°C and +40°C
- Maximum humidity at place of installation: 60%
- Limit value for dust content in ambient air: 10mg/m<sup>3</sup>

The **SST-UL** recirculation fan is used to extract and clean the exhaust air from hazardous material, safety, chemical and gas cylinder cabinets. The device ventilates the cabinet depending on the required air change rate, filters out hydrocarbons and indicates the current function by means of two lamps each (RED = Error, GREEN = Service). The concentration of hazardous substances in the exhaust air is shown in 25% steps. A concentration of 100% corresponds to the maximum load. The DSA (dynamic saturation adjustment) function corrects the switching point according to the ambient air quality. The device also has two potential-free outputs for these signals. It is suitable for extracting air from cabinets as defined in EN 14470-1 and DIN 12925 Part 2 and conforms to the requirements of TRbF 20, TRGS 510 and GefStoffV. The user is responsible for ensuring safe operation of the system in accordance with TRbF 20 Section 15. The statutory limit values (TRGS 900 "Occupational exposure limit values" and TRGS 903 "Biological limit values") must be observed at the place of use.

#### PERSONAL INJURY AND MATERIAL DAMAGE

Improper use could result in broken impellers, fatigue failures or fires due to sparks, possibly leading to personal injury or material damage.

### 12.6.2. Specific hazards

The device conforms to 2014/34/EU (ATEX) (equipment group II, category 3 and temperature classification T4); it is exclusively intended for extracting air from Zone 2 or non-Ex zones and must be installed in a safe area.

Marking: CE  II \_/3G c e IIB T4

All forms of use which deviate from or exceed the limits of use described above are considered to be contrary to the intended purpose. The manufacturer / supplier is not liable for any damages resulting from such use; the sole responsibility lies with the user.

No specific hazards for people or the environment arise provided this equipment is used as intended and only with original accessories, and provided this Installation and Operating Manual as well as the information on safe commissioning are observed.

### 12.6.3. Operating principle, safe installation, commissioning, adjustment, installation / removal, maintenance

The device is shipped ready to use and install with the exception of the control cables, which are not included in the scope of supply because they must be adapted to the specific conditions on site. You must check the electrical specifications (voltage, frequency) and compare them with the information on the rating plate prior to connecting the device to the mains.

The device may only be installed, operated and maintained by appropriately trained, instructed and authorised persons.

- Operating principle: The air containing hazardous substances is discharged from the hazardous material cabinet via a filter system comprised of a coarse filter, a fine filter and an activated carbon filter. 99.9% of all hydrocarbons are retained in this way by adsorption. This micro-controlled device monitors the volume flow and saturation of the activated carbon filter. The operational status is indicated by three LEDs (GREEN, YELLOW and RED) as well as by an audible signal. To meet the requirements of TRbF 20 Annex L, the exhaust function is interrupted as soon as the saturation LED shows "Error". This ensures that no hazardous substances are discharged from the cabinet into the workroom. The interior of the device is automatically flushed in the event of an error in order to reduce the concentration. If flushing does not cause the concentration to fall below the setpoint, the device goes into "alarm". Alarms can only be cancelled by authorised persons. Short concentration peaks (less than 9 seconds) are simply indicated (all saturation lamps on) without the device going into "error". The coarse and fine filters significantly extend the life of the activated carbon filter by retaining dust particles. The actual saturation up to the limit value is indicated by the three yellow LEDs in 25% steps. The quality of the ambient air is shown by another yellow LED. The brightness of this LED varies according to the amount of the hazardous substances.

- Installation location: the capacity of the recirculating air filter is intended for cleaning the air volume from the cabinet (cabinet volume x air turnover). Since the exchange air is taken from the ambient air via cabinet air intake, it must be taken into account that the concentration of pollutants is influenced by the composition of the ambient air. The total concentration of pollutants before the filter is therefore composed of the sum of the cabinet and ambient concentrations. When installing the cabinet, make sure that the ambient air is pollutant-neutral. The device provides an indication of the quality of the ambient air using the four-level "DAS" brightness display (in the normal state, the LED does not light up). The brightness levels are visible in test mode - operation of the test button.  
Polluted air circulation can result in frequently occurring fault messages and can eventually lead to shutdowns due to increased concentration values. In order to ensure the normal operation of the system, the following installation situations should be avoided:
  - Proximity of cleaning facilities where solvent or cleaning agents are handled.
  - Proximity of storage of outgassing substances.
  - Proximity of workplaces that generate a lot of dust (sanding, polishing, sand blasting etc.)
  - Proximity of workplaces that generate a lot of moisture (washing facilities, cooking facilities, chemical/thermal surface treatment baths, dryers, etc.)
  - Proximity of nearby painting activity
  - Proximity of exhaust gas emissions (engine tests in repair shops, forklift traffic, incinerators etc.)
- Installation and commissioning: The device is attached to the cabinet by means of the connection fitting (accessory). The potential-free signal outputs are connected by the customer if necessary. The device goes into test mode as soon as it is connected to the mains supply with the mains plug. The LEDs, potential-free outputs and sensors are then checked one at a time.
- Function check: A function check should be carried out whenever the device is started up or as needed (see "Troubleshooting"). When the device is connected, or if you press the TEST button, all lamps light up for 3 seconds. Following this, the relays for the potential-free contacts are switched on and off for 3 seconds each. The flow and saturation sensors are checked in the final test sequence. The flow is checked by switching the technical ventilation system on and off (5 seconds on, then 5 seconds off). The LEDs for the flow change colour from RED to GREEN and the result of the test is indicated during the pauses (RED = Error, GREEN = OK). At the end of the test, the device automatically goes into normal mode. The exhaust air monitoring function can also be checked by briefly covering over the exhaust vent in order to throttle the volume flow.
- Adjustment: The exhaust air monitoring function is set to the required flow and saturation at the factory, so that no further adjustments are necessary when you start up the device for the first time.

- Repairs: Repair work may only be carried out by authorised persons.
- Maintenance: The device must be checked by an authorised person as and when necessary (once a year as a minimum) in accordance with TRbF 20 Section 15.4. The filters must be replaced or cleaned at the same time. Replacement filters can be ordered from the supplier. We recommend stocking a set of replacement filters to avoid interrupting operation. The device records the operating time internally and outputs an audible signal (one 3-second beep every hour after 365 days) to indicate that maintenance is due.


**12.6.4. Troubleshooting**

Fault	Cause	Remedy
Device not working	No power Fuse defective Device defective	Check the mains voltage Replace the fuse of the mains connector Send in the device
Red flow lamp ON	Dust filter clogged  Volume flow interrupted	Replace the coarse filter Replace the fine filter Remove the obstruction from the exhaust vent
Red saturation lamp ON, audible signal, air flow interrupted	Activated carbon filter saturated Concentration of hazardous substances too high  Other hazardous substances contained in exhaust air	Replace the filter Check the contents of the cabinet and remove any leaked substances Start up the device again (authorised persons only) Contact the supplier

### 12.6.5. Specifications


Model	: SST-UL	
Ventilation system	: RCRv75/2m	
Flow	: 10 to 25 m <sup>3</sup> /h	
Total pressure	: Max. 250 Pa	
Sound level	: 38 dB(A) at a distance of 3 m	
Monitored parameters	: - Flow	<div style="border: 1px solid black; padding: 5px; width: fit-content;"> <p style="text-align: center; margin: 0;"><b>Test</b></p> <p style="text-align: center; margin: 0;">Volumenstrom Flow</p> <div style="display: flex; justify-content: space-around; margin: 0;"> <span>○</span><span>○</span><span>○</span> </div> <p style="text-align: center; margin: 0;">NC COMM NO</p> <p style="text-align: center; margin: 0;">Sättigung Saturation</p> <div style="display: flex; justify-content: space-around; margin: 0;"> <span>○</span><span>○</span><span>○</span> </div> <p style="text-align: center; margin: 0;">NC COMM NO</p> <p style="text-align: center; margin: 0;">230 V / 50 Hz</p> <p style="text-align: center; margin: 0;">In = 0,2 A</p> <p style="text-align: center; margin: 0;">Fuse 0,6 A</p> </div>
- Saturation of activated carbon filter		
LEDs	: ROT = Error, GREEN = Service	
	: 25 / 50 / 75%	
	: Ambient concentration	
Signal outputs	: - Potential-free for flow	
- Potential-free for saturation of activated carbon filter		
Contact power	: 230 V / 8 A	
Power supply	: 230 V / 50 Hz	
Current consumption	: Max. 0.5 A	
IP protection	: CE  II _/3G c e IIB T4	
Frequency	: 50 Hz	
Motor protection	: PTO 90°C	
Weight	: 20 kg	

### 12.6.6. Characteristic

**12.6.7. Dimensions**

**12.6.8. Accessories shipped with the SST-UL**

**12.6.9. Installing the SST-UL**

- Insert the adapter
- Place the device on top of the cabinet with the adapter over the exhaust vent, then place the height adjuster underneath the device and align
- Fasten the height adjuster to the cabinet using the screws provided

## 12.7. Exhaust fan, type SST AL

### 12.7.1. Instructions for safe use

The **SST-AL** exhaust fan is used to extract hazardous substances from hazardous material, safety, chemical and gas cylinder cabinets as well as from acid and lye cabinets (aggressive and / or explosive media with a dust content of  $< 5 \text{ g/m}^3$ ). Ambient temperatures from  $+10^\circ\text{C}$  to  $+40^\circ\text{C}$  are permitted. The device ventilates the cabinet depending on the required air change rate and indicates the current function by means of two lamps each (RED = Error, GREEN = Service). The device also has potential-free outputs for these signals. It is suitable for extracting air from cabinets as defined in EN 14470-1 and DIN 12925 Part 2.

#### **PERSONAL INJURY AND MATERIAL DAMAGE**

Improper use could result in broken impellers, fatigue failures or fires due to sparks, possibly leading to personal injury or material damage.

### 12.7.2 Specific hazards

The device conforms to 2014/34/EU (ATEX) (equipment group II, category 3 and temperature classification T4); it is exclusively intended for extracting air from Zone 2 or non-Ex zones and must be installed in a safe area. Marking: CE  II \_/3G c e IIB T4

All forms of use which deviate from or exceed the limits of use described above are considered to be contrary to the intended purpose. The manufacturer / supplier is not liable for any damages resulting from such use; the sole responsibility lies with the user.

No specific hazards for people or the environment arise provided this equipment is used as intended and only with original accessories, and provided this Installation and Operating Manual as well as the information on safe commissioning are observed.

### 12.7.3 Safe installation, commissioning, adjustment, installation / removal, maintenance

The device is shipped ready to use and install with the exception of the control cable and exhaust duct, which are not included in the scope of supply because they must be adapted to the specific conditions on site. You must check the electrical specifications (voltage, frequency) as well as the resistance of the system and compare them with the information on the rating plate prior to connecting the device to the mains.

The device may only be installed, operated and maintained by appropriately trained, instructed and authorised persons.

- **Installation and commissioning:** The device is attached to the cabinet by means of the connection fitting (included). The signal output and the exhaust duct are connected. You must ensure that no foreign objects can get into the exhaust air monitoring system via the exhaust duct (protective grid). The device is ready to use as soon as it is connected to the mains supply with the mains plug. The green lamp indicates normal mode.
- **Function check:** A function check should be carried out prior to connecting the exhaust duct or as needed (see "Troubleshooting"). To do this, attach the device and connect it to the mains supply. The device is ready to use as soon as it is connected. The green lamp now lights up permanently. The device should go into error mode if you cover over the exhaust air outlet briefly in order to throttle the volume flow. Repeat this procedure several times. In order to test the monitoring function without removing the device, press the TEST button on the rear.
- **Adjustment:** The exhaust air monitoring function is set to the required flow at the factory, so that no further adjustments are necessary.
- **Installation / removal:** Repair work may only be carried out by authorised persons.
- **Maintenance:** The device requires no maintenance.


### 12.7.4 Troubleshooting

Fault	Cause	Remedy
Device not working	No power Device defective	Check the mains supply Send in the device
Red lamp permanently ON or blinking, function check OK	System resistance too high	Change the duct layout, reduce the air velocity in the duct (by increasing the cross-section)
Red lamp permanently ON, function check not OK	Ventilation system sticks, control loop defective	Send in the device

### 12.7.5 Specifications


Model	: SST-AL	
Ventilation system	: RCRv75/2m	
Ventilation housing material	: Acid and lye-resistant plastic	
Impeller material	: Acid and lye-resistant plastic	
Flow	: 50 to 200 m <sup>3</sup> /h	
Total pressure	: 250 to 50 Pa	
Speed	: 2400 to 2800 rpm	
Sound level* Lp2A	: 35 dB(A)	
LEDs	: ROT = Error, GREEN = Service	
Operating cycles	: Min. 1,000,000	
Signal output	: Potential-free via DIN bush	
Contact power	: Max. 230 V / 8 A	
Power supply	: 230 V +/-10%	
Current consumption	: Max. 0.45 A	
Frequency	: 50 Hz	
Motor protection	: PTO 90°C	
IP protection	: CE  II _/3G c e IIB T4	
Weight	: 8 kg	

**Signal output assignment**  
Max. 230 V / 8 A


\* At a distance of 3 m in accordance with DIN 45635, with air connections on the suction and discharge sides

12.7.6 Characteristic

12.7.7 Dimensions

### 12.7.8 Accessories shipped with the SST-AL

	1. Operating Manual
	2. Connectors for potential-free contacts
	3. Connection fitting for all cabinet types, already fitted
	4. Connection fitting (type SST-D only)
	5. Adapter for third-party equipment
	6. Adapter for SST-P

### 12.7.9 Installing the SST-AL for cabinets

- Insert the adapter (6) at the free end of the connection fitting and fasten it with the hose clamp provided
- After preparing the adapter in this way, insert it into the device and the exhaust air outlet of the cabinet and fasten it to the cabinet outlet using the hose clamp provided. Align the device

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**Information on the GS Mark Certificate**

This GS Mark Certificate only applies to the company on the first page and the manufacturing facility stated there. It can be transferred to third applicants by the certification body only.

Only the above stated product is allowed to mark with the GS Mark as depicted on page 1.

The GS Mark Certificate can be declared invalid or rather it can be withdrawn (paragraph 21 clause 5 ProdSG).

The certification body inspects the manufacture of the product which is specified in this certificate and keeps the lawful use of the GS Mark under surveillance.

The manufacturer is obliged to meet the conditions which are necessary for the correct manufacture. For it has among others continually to monitor the manufacture of the products fitted out with the GS Mark on conformity with the examination specifications. Further more he has to examine or rather to tolerate the monitoring examinations which laid down in the examination specifications or required by the certification body.

The test and certification directives and the terms and conditions of the acknowledged body apply to this GS Certificate.

**Protectoplus**

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Grüner Kamp 19-21

D-24768 Rendsburg

Telefon +49 4331-45 16 0

Telefax +49 4331-45 16 11

[info@protecto.de](mailto:info@protecto.de)

[www.protecto.de](http://www.protecto.de)

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MPA Dresden GmbH Geschäftsführer: Thomas Hübler Amtsgericht Chemnitz HRB 28268 Sparkasse Mittelhachsen