CEMO Sicherheitsschrank

D
2-19

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(Unterschrift)

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Operating instructions

- provide to operator.
- must be read before using the equipment for the first time.
- retain for future use.

Important!

These operating instructions must be stored near the safety cabinet in an environment protected from moisture and heat.

Important!

Do not amend the contents of these operating instructions. Do not damage, amend or remove any part of these operating instructions.

If required, new operating instructions for the safety cabinet can be obtained from manufacturer CEMO.

Important!

These operating instructions must be passed on with the safety cabinet if it is sold.

Dear Customer,

Thank you for choosing a quality product from CEMO.

Our products are manufactured using modern production methods and are subject to quality control measures. We do everything we can to ensure that you are satisfied with our product and that it is user friendly.

If you have any questions about your product, please get in touch with your dealer or contact our Sales department directly.

Kind regards,

Eberhard Manz, Managing Director

| | | | |
|--|-----------|---|-----------|
| Operating instructions | 20 | | |
| 1. General provisions | 22 | 15. Accessories | 32 |
| 1.1 Safety | 22 | 15.1 Spare keys | 32 |
| 1.1.1 Maintenance and monitoring | 22 | 15.2 Bottom collecting trays with perforated plate shelf | 32 |
| 1.1.2 Using original parts | 22 | 15.3 Tray bottom | 32 |
| 1.1.3 Operating the safety cabinet | 23 | 15.4 Radial fans | 32 |
| 1.2 Intended use | 23 | 15.4.1 Radial fan model 1 | 32 |
| 1.2.1 Summary | 23 | 15.4.2 Radial fan model 2 | 32 |
| 1.3 Inappropriate use | 23 | 15.4.3 Radial fan model 3 | 32 |
| 2. Safety instructions | 23 | 15.4.4 Accessories for radial fans: | 32 |
| 2.1 Safety precautions | 23 | 15.5 Exhaust air fan with airflow monitoring | 33 |
| 2.2 Type plate | 23 | 15.6 Recirculation fan with activated carbon filter and airflow monitoring | 33 |
| 2.3 Disclaimer | 23 | 15.7 Cable grommet for safety cabinet PROline 6/20 & 12/20 | 33 |
| 2.4 General safety instructions | 23 | 15.8 "Active storage" accessories package | 33 |
| 2.5 Instructions for storage | 24 | 16. Disposal | 34 |
| 2.6 Legal requirements | 24 | 17. Warranty | 34 |
| 2.7 Product and manufacturer identification | 24 | 18. Drawings | 35 |
| 3. Models | 25 | 19. Technical data | 36 |
| 4. Technical specifications | 25 | 20. Declaration of conformity | 37 |
| 5. Transport | 25 | | |
| 5.1 Dismantling the transport packaging | 25 | | |
| 6. Installation | 25 | | |
| 7. Commissioning | 26 | | |
| 7.1 Installing attachments | 26 | | |
| 7.1.1 Ventilation nozzles | 26 | | |
| 7.1.2 Cable grommet (optional) | 26 | | |
| 7.1.3 Fans (optional) | 26 | | |
| 7.2 Aligning the cabinet | 26 | | |
| 8. Interior fittings | 27 | | |
| 8.1 Tray bottoms (height-adjustable) | 27 | | |
| 8.4 Door closure | 27 | | |
| 9. Storage quantities | 28 | | |
| 11. Grommets | 29 | | |
| 12. Earthing | 29 | | |
| 13. Decommissioning after a fire | 30 | | |
| 14. Maintenance and servicing | 31 | | |
| 14.1 General | 31 | | |
| 14.2 Annual safety inspection | 31 | | |
| 14.3 Cleaning | 31 | | |

1. General provisions

The safety cabinet is state-of-the-art technology and complies with recognised technical safety regulations.

The safety cabinet carries the CE symbol, i.e. the EU directives and harmonised standards of relevance to the safety cabinet were applied during design and manufacture. The fire-resistant construction has been successfully tested in accordance with EN 14470-1.

The safety cabinet must only be used in a technically flawless condition in the version delivered by the manufacturer. For safety reasons, do not make unauthorised modifications to the safety cabinet.

1.1 Safety

Each safety cabinet is tested for functionality and safety before delivery.

When used as intended, the safety cabinet is safe to operate.

Incorrect operation, failure to observe the safety instructions or misuse may result in risks to:

- the life and limb of the operator
- the safety cabinet and other assets of the operator.
- the function of the safety cabinet.

As the operator of the safety cabinet, you are responsible for ensuring that

- all safety instructions are understood and observed.
- compliance with the applicable work safety and fire protection regulations.
- only instructed individuals operate the safety cabinet.
- the opening area of the doors is kept clear.
- the doors are always kept closed.
- the doors are locked with the keys provided to prevent unauthorised access.
- any hazardous substances that are leaking out are immediately picked up and disposed of properly.

This cabinet is not suitable for the storage of containers containing aggressive chemicals (acids/alkalis).

1.1.1 Maintenance and monitoring

The maintenance schedules proposed in these operating instructions represent the absolute minimum needed to ensure the safety and service life of the cabinet under normal operating conditions. Watch out for any malfunctions or potential safety issues at all times during operation. In the event of a malfunction, the safety cabinet must be taken out of operation immediately!

Obligations of the user:

- Compile the operating instructions and publish them in an appropriate place in the workplace.
- Carry out a risk assessment.
- Identify activities of named personnel.

The safety cabinet must be regularly checked to ensure it is in a safe condition. (see section 13).

This check shall include:

- A visual inspection for damage (seals and door closing mechanism, etc.)
- Functional testing
- Completeness/identifiability of the warning, mandatory and prohibitory signs on the safety cabinet
- The prescribed maintenance work.

1.1.2 Using original parts

Only use original parts provided or recommended by the manufacturer. Take note of all safety and usage information provided with these parts.

This concerns spare and wear parts. For safety reasons, making any modifications to the safety cabinet is not permitted.

Important!

Any unauthorised modification to this product without the written consent of the manufacturer, CEMO, will invalidate the warranty.

1.1.3 Operating the safety cabinet

The safety cabinet must only be operated by instructed persons who

- have read and understood the operating instructions.
- have proven their ability to operate the equipment.
- have been assigned to operate the equipment.

Important!

These operating instructions must be available at the safety cabinet so that they are easily accessible to all users.

1.2 Intended use

The safety cabinet is suitable for the storage of flammable liquids in closed containers at the workplace (in work rooms) in compliance with national rules and regulations.

1.2.1 Summary

Any other use is considered inappropriate. The intended use includes compliance with all the information in these operating instructions.

1.3 Inappropriate use

Important!

Inappropriate use also includes failure to comply with the directions in these operating instructions.

It also includes:

- Non-compliance with the provisions of national law that apply in each case.
- Storage of media other than those specified under intended use.

2. Safety instructions

2.1 Safety precautions

Inappropriate use or installation of this product may result in serious injury or death.

- Read and follow all warnings and precautions for safe use.
- Servicing, maintenance and inspections must be carried out by qualified personnel.

2.2 Type plate

A type plate containing the essential product data is attached in a clearly visible position on the front side of the safety cabinet. This type plate must not be removed.

2.3 Disclaimer

Any use beyond the intended use constitutes misuse. CEMO is not liable for any damage resulting from misuse.

CEMO has no influence over the liquids used by the operator. The user must check that the cabinet is suitable for the intended application.

CEMO is not liable for damage of any kind caused by the goods stored within it.

2.4 General safety instructions

Follow the safety instructions in these operating instructions in order to reduce risks to health and avoid hazardous situations. Not using the device as intended in accordance with these operating instructions poses the risk of accidents and inadequate fire protection.

Carry out a risk assessment:

- Select a practical location to install the cabinet (see section 6).
- Only store flammable liquids with a hazard assessment.
- Lock the cabinet to prevent unauthorised access. Always keep the doors closed. Always keep the opening area of the doors clear.

The doors are equipped with a self-closing door locking system. This is triggered at a temperature from 40°C – 50°C. The closing force, closing speed and door lock must only be changed by an authorised service technician!

2.5 Instructions for storage

Observe the national regulations, rules and laws in place for handling flammable liquids.

According to DIN EN 14470-1, the bottom collecting tray must have a collection volume of 10% of the total containers stored or 110% of the largest container (depending on which volume is greater). Therefore, do not use larger or too many containers. Observe the maximum storage quantities! Any hazardous substances that are leaking out must be picked up immediately and disposed of properly.

General provisions:

- Always keep the doors closed.
- Close drawers after use.
- Do not store any other flammable objects (packaging, etc.) in the cabinet.

In the event of a fire:

- Do not open the cabinet doors (risk of explosion and spread of fire).
- Alert the relevant authorities (for example the fire service).
- Leave the building immediately.

Please note:

Wear personal protective equipment.

When handling hazardous liquids, always wear protective equipment suitable for the work being carried out. Always wear the correct protective equipment during installation, use and maintenance: Gloves, safety glasses, safety shoes and long work clothes

2.6 Legal requirements

According to the German Occupational Health and Safety Act (ArbSchG), any hazards that can arise from flammable liquids must be assessed or evaluated in a hazard assessment. Pay particular attention to the latest version of the following regulations:

- The requirements of the fire authority
- General building law and building requirements
- Explosion protection guidelines
- Ordinance on Hazardous Substances (GefStoffV)
- Technical Rules for Hazardous Substances (TRGS, in particular TRGS 510)
- Technical Rules for Operational Safety (TRBS)
- Accident prevention regulations
- German guideline for workplaces, ArbSchG
- German Social Accident Insurance (DGUV) information
- DIN EN 14470-1

The applicable national regulations must also be observed.

2.7 Product and manufacturer identification

Name and address of the manufacturer:

CEMO GmbH

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D-71384 Weinstadt, Germany

Phone +49 7151 9636-0

Fax +49 7151 9636-98

www.cemo.de

kontakt@cemo.de

Product identification: CEMO PROline safety cabinet for flammable liquids according to DIN EN 14470-1.

3. Models

Type 90: PROline 12/20 and 6/20.


Standard equipment: Three (3) height-adjustable tray bottoms, one (1) collecting tray, one (1) perforated plate cover as a first storage level and one removable base cover.

4. Technical specifications

See 18. Drawings on page 35 and 19. Technical data on page 36.

5. Transport

- For safety reasons, the cabinet must only be transported when it is empty and closed.
- Please be aware that the cabinet is very heavy.
- There is a risk of crushing. Wear safety shoes when transporting/moving the cabinet.
- Use technical aids for transport, for example pallet trucks and forklift trucks. With the base cover removed, slide the forks under the cabinet. Also ensure that the cabinet is transported without any risk of slipping.
- When lying down, the safety cabinet can only be transported on its side or rear wall.
 - ▶ The cabinet must never be transported with the doors facing downwards!

The safety cabinet should be transported upright. Tilting must be avoided during pick-up and transport. Tilting can result in damage, which restricts or destroys the functions of the built-in safety technology. Furthermore, it is essential to avoid setting down the cabinet abruptly. Safety cabinet 6/20 can only be picked up using a fork on a forklift. There is a risk of tipping if handled improperly.

5.1 Dismantling the transport packaging


Transport locks are installed to protect the cabinet from damage during transport. The transport locks in the door gaps must only be removed once at the final installation location!

6. Installation

Important!

The cabinet may only be installed inside a building.

- Do not install it near emergency escape routes.
- Ideally, choose a completely flat space in the building with a load-bearing surface for installation (be aware of the heavy own weight of the cabinet!). The base height adjusters can be used to compensate for a difference of up to 10 mm between the highest and lowest base height adjusters.
- The operating temperature is between 5°C and 35°C.
- The installation location should be in a well-ventilated room. The cabinet must be protected against moisture.
- Fire risk. For safety reasons, the top of the cabinet should be kept clear. Exception: Fan including accessories for cabinet ventilation.

7. Commissioning

Before initial commissioning, the safety cabinet must be inspected for any damage or defects.


Important!

Refilling or filling operations (active storage of hazardous substances) must only be carried out with an additional earthed interior (see "Active storage" accessories package).

7.1 Installing attachments

7.1.1 Ventilation nozzles

The two ventilation nozzles must be attached to the ceiling of the cabinet using the screws provided. The ventilation nozzles are components relevant for fire protection and close the exhaust and supply air openings in the event of a fire at a temperature of 70°C (±10°C).

Important!

Operation without the ventilation nozzles screwed on is not permitted!

7.1.2 Cable grommet (optional)




For installation of the cable grommet, see the instructions in the accessories package.

7.1.3 Fans (optional)


For installation of the fans, see the instructions in the accessories package.

7.2 Aligning the cabinet


To ensure proper use, the cabinet must be placed flat/horizontal.

The cabinet can be leveled by adjusting the base height adjusters on the bottom of the cabinet using a size 5 Allen key. A maximum of three base height adjusters should be adjusted for this purpose. The maximum adjustment height is 10 mm.

Then check the alignment of the cabinet doors. If necessary, correct using the base height adjusters. Finally, mount the base cover on the cabinet.

The installation of an anti-tilt/wall attachment may be necessary depending on local circumstances and installation conditions (see operator's risk assessment).

8. Interior fittings



8.1 Tray bottoms (height-adjustable)

Important!

The top storage level must not be more than 1.75 m above the floor!


Installation of the tray bottoms:

The height of the tray bottoms can be adjusted (in 32 mm intervals):
Observe the permissible load of the tray bottoms (see technical data)!

8.3 Bottom collecting tray

Please note that the sealing profile of the bottom tray is always in contact with the side walls and the rear wall and is firmly seated on the edges of the bottom tray. This ensures that all liquids can be collected using the tray in the event of a leak. If the perforated plate cover is inserted into the bottom tray, it can also be used as a storage space:


Important!

The bottom tray must not be used as a storage space without a perforated plate cover!

8.4 Door closure


The safety cabinet is equipped with a door locking system, which automatically closes the doors in the event of a fire.

The doors can be locked in the open position. In the event of a fire, a built-in thermocouple ensures that the lock is released at a temperature of 40°C – 50°C and that the doors are closed as a result.

Important!

The working area of the door locking system must always be kept clear and must not be blocked by containers. The receptacles for the door locking bolts must be kept free of contamination at all times.


The hinged doors of the safety cabinet lock automatically when closed. To open the doors, the closing mechanism must be actuated using the key.

9. Storage quantities

The national regulations must be observed.

"The bottom collecting tray must have a minimum collection volume of 10% of all vessels stored in the cabinet, or at least 110% of the volume of the largest single container, whichever is the greater." (EN 14470-1)

See the Technical Rules for Hazardous Substances:
TRGS 510.

10. Ventilation

The safety cabinet has a supply air and an exhaust air connection. In the event of a fire, the supply and exhaust air openings of the safety cabinet close automatically at an airflow temperature of 70°C ($\pm 10^{\circ}\text{C}$).

Please note:

In general, connecting the safety cabinet to a ventilation system is recommended.

The position of the supply and exhaust air connections is shown in Annex 1. The exhaust air vent (DN75) must be connected to the exhaust pipe (suction side). The exhaust air must be routed via a pipe to a non-hazardous location outside the building. The supply air is extracted from the surrounding area (installation space in the building). Proper connection to a ventilation system must be checked after installation using a smoke tube, for example. The national rules, laws and regulations must be observed with regards to air changes and exhaust air routing.

In general, when the door is closed, an air change corresponding to at least 10 times the empty volume of the cabinet interior per hour must take place. The pressure difference between the pressure at the connection point of the ventilation system and the ambient pressure must not exceed 150 Pa (also see Annex 2).

Observe the position of the fans on the cabinet ceiling (figure) and the maximum load on the cabinet ceiling (20 kg), install the fan and select an appropriate fan (see accessories).

It is also possible to use the cabinet without mechanical exhaust air. However, it is important to note the following:

- The cabinet must be connected to earth to prevent electrostatic charging (see section on earthing).
- There is a zone 1 explosive area inside the cabinet.
- There is a zone 2 explosive area within a radius of 2.5 m around the cabinet and up to 0.5 m above the floor. This area must be marked in accordance with the Accident Prevention Regulations (BGV A8):

- In general, the national rules, laws and regulations must also be observed here.

Important!



A mechanical ventilation system is mandatory for the storage of highly flammable liquids (H224).

- In mechanically ventilated work rooms with at least five air changes per hour, the zone 2 explosive area can be reduced to 1 m in front of the safety cabinet and 0.5 m at the side of the safety cabinet and to a height of 0.3 m above the floor.

11. Grommets

Cable grommets can be installed in the cabinet ceiling; see accessories for the order number. The positions are marked on the cabinet as ①.

| Model | Number of grommets |
|---------------|--------------------|
| 12/20 Type 90 | 2 |
| 6/20 Type 90 | 1 |

12. Earthing

The safety cabinet has a standard earth connection ② on the ceiling.

Equipotential bonding must be installed by qualified personnel. Check the conductivity after installation!

13. Decommissioning after a fire

In the event of a fire:

- Keep calm, leave the building and call the fire department immediately.
- A flammable gas-air mixture may have formed inside the cabinet as a result of the fire. Remove all ignition sources within 10 metres and only use non-sparking tools.
- The cabinet must only be opened by qualified personnel (the fire service) after a minimum of 24 hours have passed. Take extreme care when opening the doors after a fire. Wear personal protective equipment and have suitable extinguishing agents ready.

Important!

Do not open the cabinet if the cabinet surface is still warm. The cabinet should only be opened outside and by qualified personnel (for example the fire service). Opening the doors of the cabinet too early within a building can cause fires to spread and put lives at risk. In the event of a fire, ensure that personal protective equipment (PPE) is worn.

Do not reuse the cabinet if it has been damaged as a result of a fire or the use of extinguishing agents.

14. Maintenance and servicing

14.1 General

The cabinet is considered to be a safety system (in accordance with Section 4, paragraph 3 of the German Workplace Ordinance, c, EEC Directive 89/391 and Rule No. 108-007 of the German Social Accident Insurance Association) and must therefore be checked regularly for safety and functionality. The safety cabinet must always be checked for externally visible faults or damage in the following cases:

- Prior to use,
- Following changes,
- After maintenance work.

If damage or faults are visible, the cabinet must be taken out of service until these faults have been rectified. The following maintenance work must be carried out at certain intervals:

Important!

Spilled liquids must be removed immediately and disposed of properly.

| Interval | Assembly | Action |
|-------------|-----------------|--|
| As required | Doors | Lubricate the door latch and hinges. |
| Daily | Collecting tray | Ideally daily, but at least once a week, the collecting tray should be inspected for any leakage of liquids. Any liquid that has leaked out must be removed immediately and the collecting tray and the platforms must be inspected to establish whether it is still suitable for use; it must be replaced if necessary. |
| Monthly | Cabinet | Clean and check for damage. |
| | Door closure | Open the door and check that it is closing properly. |
| | Ventilation | Test the effectiveness using a wool thread, meter, smoke tube, etc. |
| | Seals | Check for correct positioning and replace the seals if they are visibly damaged |
| Annually | Cabinet | Inspection by a specialist |

For collecting trays compliant with StawaR (Steel Tray Directive): "The condition of the collecting tray in safety cabinets according to DIN EN 14470-1 must be checked every 6 months by way of a visual inspection. At least in the case of steel collecting trays in accordance with DIN EN 10025-2 or DIN EN 10028-2, the test must also be carried out on the underside of the tray. The result must be recorded and presented to the responsible water authorities upon request."

14.2 Annual safety inspection

According to Section 4, paragraph 3 of the German Workplace Ordinance (ArbStättV) and BGR 234/6.1 of the German Employers' Liability Insurance Association Rules for Safety and Health at Work, the safety and functionality of safety cabinets must be checked.

The outcomes and measures must be documented.

The annual inspection can only be carried out by a specialist. Only CEMO can guarantee optimum servicing of CEMO products.

We are happy to carry out the annual inspection of your safety cabinet, including as part of a maintenance contract.

Send inquiries to:

service@cemo.de or call

Tel: +49 7950 /9803-2222

The due date for the next annual inspection can be found on the inspection sticker on the type plate affixed to the cabinet.

14.3 Cleaning

The safety cabinet can be cleaned using a soft cloth.

15. Accessories

15.1 Spare keys

Spare keys for door lock, simultaneous locking, order number on request

15.2 Bottom collecting trays with perforated plate shelf

Bottom collecting tray, 66 litres with perforated plate shelf for safety cabinet PROline 12/20 (order no. on request)

15.3 Tray bottom

Additional tray bottoms can be retrofitted:

Tray bottom for safety cabinet PROline 12/20, order no. 11881

Tray bottom for safety cabinet PROline 6/20, order no. 11912

15.4 Radial fans

Radial fans must be installed by specialist personnel (note: not ready to plug in!)

15.4.1 Radial fan model 1

Order no. 11431

ATEX-compliant:

Extraction from zone 2 and installation in a non-hazardous area. CE Ex II 2G Ex h IIB T3 Gc

15.4.2 Radial fan model 2

Order no. 11432

ATEX-compliant:

Extraction from zone 1 and installation in zone 1.

Motor: II 2G Ex eb IIC T4 Gb, fan: Ex II 2G Ex h IIB T4 Gb

15.4.3 Radial fan model 3

Order no. 11433

ATEX-compliant:

Extraction from zone 1 and installation in

zone 1. Motor: II 2G Ex eb IIC T4 Gb, fan: Ex II 2G Ex h IIB T4 Gb

15.4.4 Accessories for radial fans:

| Designation | Order no. |
|---|-----------|
| Throttle valve, DN 75, PPs | 11434 |
| Connection nozzles, NW 120-75, including clamp | 11435 |
| Extension piece, sleeve, DN 75 | 11436 |
| Extension piece, connector, DN 75 | 11437 |
| Ventilation hose L=750 mm, DN 75, including 2 clamps 11438 | 11438 |
| Ventilation hose L=2000 mm, DN 75, including 2 clamps 11439 | 11439 |
| T-piece, DN 75 11440 | 11440 |

15.5 Exhaust air fan with airflow monitoring

Order no. 8739

Solution ready to plug in.

For extraction from explosion protection zone 2

and installation in explosion-free zone

Protection class: CE Ex II _/3 G IIB T4

15.6 Recirculation fan with activated carbon filter and airflow monitoring

Order no. 10963

Solution ready to plug in.

ATEX-compliant:

For extraction from explosion protection zone 2

and installation in explosion-free zone

Protection class: CE Ex II _/3 G IIB T4

15.7 Cable grommet for safety cabinet PROline 6/20 & 12/20

In the ceiling of the cabinet, for self-assembly, sheathed cables, max. 5 x 2.5 mm² cross-section
Order no. 11916

15.8 "Active storage" accessories package

During active storage (refilling, filling of hazardous materials in the safety cabinet), the interior of the cabinet must be earthed.

Order no. on request.

16. Disposal

The safety cabinet can be dismantled into individual parts and recycled.

Pure gypsum fibreboards and plasterboards are considered to be construction and demolition waste and do not contain any hazardous substances.

All plastic and non-degradable material parts must be collected separately and recycled by an authorised disposal firm.

Dispose of scrap metal at a scrap metal recycling centre. In general, the national and local disposal regulations must be observed.

17. Warranty

We guarantee that the safety cabinet will be produced free from defects in functionality or workmanship under our general terms and conditions of trade.

These can be viewed at
<https://shop.cemo.de/agb/>


The warranty applies only under the condition that the above operating and maintenance instructions and all applicable regulations are closely followed. Modification of the safety cabinet by the customer without consultation with manufacturer CEMO GmbH invalidates any claims under the statutory warranty.

Disclaimer:

- The company "CEMO GmbH" is not liable for damage caused by inappropriate use.
- CEMO has no influence over the liquids used by the operator. The user must check that the cabinet is suitable for the intended application.
- CEMO is not liable for damage of any kind caused by the goods stored within it.

18. Drawings

Type 90: PROline 12/20 and 6/20.

19. Technical data

| Model | 12/20 Type 90 | 6/20 Type 90 |
|--|-----------------------------|---|
| Order number | 11880 | 11911 - Door hinge left 11916 - Door hinge right |
| External dimensions (WxDxH) | 120 x 60 x 196 cm | 60 x 60 x 196 cm |
| Internal dimensions (WxDxH) | 105 x 49 x 157 cm | 45 x 49 x 149 cm |
| Tare weight with tray bottoms & bottom collecting tray with perforated plate shelf | approx. 450 kg | approx. 265 kg |
| Max. total weight when fully loaded | Tare weight + 300 kg | Tare weight + 300 kg |
| Usable interior volume of the cabinet | approx. 0.81 m ³ | approx. 0.33 m ³ |
| Total storage volume in the cabinet | 300 litres | 300 litres |
| Volume of the bottom collecting tray | 33 litres | 33 litres |
| Max. volume of the largest container | 30 litres | 30 litres |
| Max. load capacity of each tray bottom (evenly distributed load) | 75 kg | 75 kg |
| Maximum load on cabinet ceiling | 20 kg | 20 kg |

20. Declaration of conformity

EC declaration of conformity

The manufacturer / distributor
CEMO GmbH
In den Backenländern 5
D-71384 Weinstadt, Germany

hereby declares that the following product

| | |
|----------------------|--|
| Product designation: | safety cabinet |
| Model designation: | CEMO |
| Type designation: | safety cabinet PROline Typ 90 6/20 and 12/20 |
| Serial numbers: | 11880, 11911, 11916 |

Description:

Safety cabinet for the storage of flammable liquids in closed containers at the workplace (in workrooms) in compliance with national rules and regulations.

complies with all relevant specifications of the applicable regulations (below), including any amendments applicable at the time of the declaration. The manufacturer bears sole responsibility for issuing this declaration of conformity.

The following harmonised standards have been applied:

EN 14470-1:2023-09 Fire safety storage cabinets - Part 1: Safety storage cabinets for flammable liquids

The following legislation has been applied:

Directive 2006/42/EC

Name and address of legal entity authorised to compile the technical documentation:

CEMO GmbH
In den Backenländern 5
D-71384 Weinstadt, Germany
Location: D-71384 Weinstadt, Germany
Date: 31/10/2023

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